


MEASURING YOUR HORSE ACCURATELY FOR REGISTRATION

HEAD LENGTH:

Place the end of a cloth measuring tape between the horse's ears at the center of the poll, allowing the tape to lay flat, down along the center of the face to an imaginary line directly above the top of the nostrils as illustrated.



CIRCUMFERENCE OF MUZZLE:

Center the end of your tape directly above the top of the nostrils, wrapping the tape around the entire muzzle and back to the starting point.



EAR LENGTH:

Measure from the tip of the ear down to the base, on the inside center of the ear as shown.




WIDTH BETWEEN EYES:

Measure across the forehead from the inside edge of one eye to the inside edge of the other eye.




LENGTH OF NECK (topline):

Begin at the center of the poll, bringing the tape down the crest of the neck to the center of the wither.



LENGTH OF BACK:

Begin at the center of the wither, following the natural contour of the back in a straight line to the beginning of the croup, at the center of the back between the points of the hip.



LENGTH OF CANNON:

Again, from the front of the leg, begin at the center of the knee and bring the tape down the lower part of the leg to the top of the fetlock joint.



CIRCUMFERENCE OF CANNON:

Measure from a point halfway between the knee and the fetlock joint, around the lower part of the foreleg.



LENGTH OF PASTERN:

Follow the line of the pastern with your tape in front, from the base of the fetlock joint to the coronet. Please note this method of measurement has been altered on the enclosed application form.



DISTANCE FROM GROUND TO BODY:

Measured from the bottomline of the heart girth (in a line directly below the base of the wither) down to the ground.

